



70% Isopropyl Alcohol Spray

Effective bactericidal action against *e.coli*, *pseudomonas aeruginosa* and a diverse range of micro-organisms.

Overview

Virox 70% Isopropyl Alcohol Spray enables an instant and efficient disinfection to a variety of hard surfaces such as tables, countertops, wheelchairs, surgery tables and x-ray equipment. Ideal for all healthcare professionals who require an effective solution to their infection control needs. Virox 70% Isopropyl Alcohol Spray is designed for hard surface disinfecting and cleaning and is meant for everyday use as well as non-smearing so they do not streak hard surfaces.

Features and Benefits



Kills 99.99% of SARS-CoV-2 in 1 minute



Designed for hard surface disinfecting and cleaning



This range carries CE Mark



Manufactured under ISO 9001, 22716 & 13485 quality system to ensure complete quality control



Effective in killing SARS-CoV-2, Influenza A H1N1 Virus



Formulated in Singapore

Active Ingredients

Isopropyl Alcohol 70% by Volume

Antimicrobial Effectiveness Tested

Test Microorganism	EN Method	Contact Time (Dirty Condition)
MRSA Escherichia coli Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae	EN13727	1 min
Candida albicans	EN13624	1 min
Virucidal against enveloped viruses (Vaccinia Virus)	EN14476	1 min
SARS-CoV-2 Influenza A H1N1	EN14476	1 min
Mycobacterium terrae	EN14348	2 min

3 products in this range

Product Code	Volume	Bottle Spray Dimensions	EAN Barcode	Units / Carton	Carton Size & Weight	Pallet Weight	Pallet Height	Cartons / Pallet	Cartons / Layer	Layers / Pallet
VR4040	100ml	D38mm x H148mm	8886389740407	12	L115 x W150 x H155 1.6kg	1045	1700	640	64	10
VR4036	500ml	D27mm x H272mm	8886389740360	12	L335 x W200 x H290 6.6kg	615	1600	90	18	5
VR4000	1L	D93.5mm x H291mm	8886389740001	6	L310 x W190 x H310 6.3kg	588	1700	90	18	5

